

AVANOS

Digestive Health
ENFit® Product Codes

CORFLO*

YOUR PREFERRED
ENTERAL PARTNER

TABLE OF CONTENTS

ENTERAL FEEDING TUBES4

CORTRAK® Enteral Access System (EAS).....	4
CORTRAK® Enteral Access System (EAS) Accessory Item	4
CORTRAK® CORFLO® Nasoenteric Feeding Tube, polyurethane, unweighted with electromagnetic transmitting stylet	4
CORFLO® Nasoenteric Feeding Tube, polyurethane, fully radiopaque, cm markings, unweighted, stylet.....	5
CORFLO® Nasoenteric Feeding Tube, polyurethane, fully radiopaque, cm markings, unweighted, non stylet.....	5
CORFLO® Nasoenteric Feeding Tube, polyurethane, fully radiopaque, cm markings, weighted, stylet.....	5
CORFLO® Nasointestinal Endoscopically Placed Feeding Tubes	6
CORGRIP® NG/NI Tube Retention System.....	6
CORFLO® PEG Kit.....	7
CORFLO® PEG Jejunal Tube, polyurethane, radiopaque, non-sterile.....	8
CORFLO® PEG Adapter Repair Kit.....	8

ENTERAL ACCESSORIES9

CORFLO® PEG Accessories	9
CORFLO® Tube Accessory	9
FARRELL® Valve System	9
aspHirate pH Indicator Strips.....	10
Enteral Syringe, single use.....	10
CORSTOP® a.c.e. stopper	11
CLOG ZAPPER® Enteral Feeding Tube Declogging System	11

ENTERAL FEEDING TUBES

CONFIRMATION WITHOUT DELAY



CORTRAK* Enteral Access System (EAS) ^Δ

PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORTRAK* Enteral Access System	20-0950	1 unit
CORTRAK* Accessory Kit	20-0951	1



CORTRAK* Enteral Access System (EAS) Accessory Item ^Δ

PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
QLN220 Printer	20-0928	1
QLN220 Printer charger	20-0929	1
Printer Paper Rolls	20-0911	3
Plastic Covers (sleeves) for SRU	20-0930	100
Stand	20-0920	1
Receiver Unit	20-0952	1
Interconnect Cable	20-0966	1

The exclusive AVANOS* CORFLO* yellow tube design offers a variety of features including:

- Unique scooped out distal tip, for easier aspiration, anti-clog
- Fully radiopaque tube with increased barium tip, for easy visualisation on x-ray
- Centimeter markings along the length of the tube
- Water-activated lubricant on the external tube tip to ease insertion
- Can stay in place as long as functional

Broad portfolio - polyurethane, longer-lasting

CORTRAK* CORFLO* Nasoenteric Feeding Tube, polyurethane, unweighted with electromagnetic transmitting stylet ^Δ



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORTRAK* CORFLO* Nasoenteric Feeding Tube 8Fr x 91cm, polyurethane, unweighted with electromagnetic transmitting stylet	42-9368TRAK2	Box of 10
CORTRAK* CORFLO* Nasoenteric Feeding Tube 8Fr x 109cm, polyurethane, unweighted with electromagnetic transmitting stylet	42-9438TRAK2	Box of 10
CORTRAK* CORFLO* Nasoenteric Feeding Tube 8Fr x 140cm, polyurethane, unweighted with electromagnetic transmitting stylet	42-9558TRAK2	Box of 10
CORTRAK* CORFLO* Nasoenteric Feeding Tube 10Fr x 91cm, polyurethane, unweighted with electromagnetic transmitting stylet	42-9361TRAK2	Box of 10
CORTRAK* CORFLO* Nasoenteric Feeding Tube 10Fr x 109cm, polyurethane, unweighted with electromagnetic transmitting stylet	42-9431TRAK2	Box of 10
CORTRAK* CORFLO* Nasoenteric Feeding Tube 10Fr x 140cm, polyurethane, unweighted with electromagnetic transmitting stylet	42-9551TRAK2	Box of 10
CORTRAK* CORFLO* Nasoenteric Feeding Tube 12Fr x 91cm, polyurethane, unweighted with electromagnetic transmitting stylet	42-9362TRAK2	Box of 10

^ΔLimited market availability, contact your sales representative

CORFLO* Nasoenteric Feeding Tube, polyurethane, fully radiopaque, cm markings, unweighted, stylet



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 5Fr x 56cm	42-9225	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 6Fr x 38cm	42-9156	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 6Fr x 56cm	42-9226	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 6Fr x 91cm	42-9366	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 8Fr x 56cm	42-9228	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 8Fr x 91cm	42-9368	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 8Fr x 109cm	42-9438	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 8Fr x 140cm	42-9558	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 10Fr x 91cm	42-9361	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 12Fr x 91cm	42-1362	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 10Fr x 109cm	42-9431	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 10Fr x 140cm	42-9551	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 12Fr x 109cm	42-9432	Box of 10

Δavailable on special order

CORFLO* Nasoenteric Feeding Tube, polyurethane, fully radiopaque, cm markings unweighted, non stylet



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO* Nasogastric/Nasointestinal Feeding Tube, 12Fr x 109cm	42-1432	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube, Bullet Bolus Tip, 6Fr x 56cm	42-1226	Box of 10

CORFLO* Nasoenteric Feeding Tube, polyurethane, fully radiopaque, cm markings, weighted, stylet



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 6Fr x 56cm	42-8226	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 6Fr x 91cm	42-7366	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 8Fr x 91cm	42-7368	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 8Fr x 109cm	42-2438	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 8Fr x 140cm	42-7558	Box of 10
CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet, 10Fr x 91cm	42-7361	Box of 10

CORFLO* Nasointestinal Endoscopically Placed Feeding Tubes

For endoscopic or radiological Over The Wire (OTW) placement;
maintains same features as standard CORFLO* Feeding Tubes

152cm CORFLO* Feeding Tube, fully radiopaque, polyurethane, regular cm, 300cm teflon coated flexible guidewire (0.035mm diameter) with dispenser sheath, and 10 Fr, 56cm length Nasal/Oral Transfer Tube



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO* Nasointestinal tube, 12Fr, 152cm, 7gr weighted	51-4602	1
CORFLO* Nasointestinal tube, 10Fr, 152cm, unweighted	51-9601	1

CORGRIP* NG/NI Tube Retention System

The CORGRIP* Nasogastric/Nasointestinal (NG/NI) Tube Retention System is indicated for use with enteral feeding tubes of 8 Fr and greater and NG decompression, suction and drainage tubes up to 18 Fr to prevent inadvertent removal or displacement of the tubes for adult patients

- Unique delivery system involves no dragging of catheter behind vomer bone
- Centimetre markings on catheters for ease and parallel insertion²
- Prevents multiple tube placements and associated discomfort and potential trauma¹
- Helps avoid the costs and risks associated with displaced feeding tubes¹
- Available for tubes of 8 FR up to 18 FR



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORGRIP* NG/NI Tube Retention System for use with 8Fr NG/NI feeding tubes	25-008	Box of 10
CORGRIP* NG/NI Tube Retention System for use with 10Fr NG/NI Feeding tubes	25-010	Box of 10
CORGRIP* NG/NI Tube Retention System for use with 12Fr NG/NI feeding tubes	25-012	Box of 10
CORGRIP* NG/NI Tube Retention System for use with 14Fr NG/NI feeding tubes	25-014	Box of 10
CORGRIP* NG/NI Tube Retention System for use with 16Fr NG/NI feeding tubes	25-016	Box of 10
CORGRIP* NG/NI Tube Retention System for use with 18Fr NG/NI feeding tubes	25-018	Box of 10



PRODUCT FEATURES

- Unique deployment - no dragging of tape around vomer bone
- Unique Slot Slide Lock (SSL) securement clip
- Each catheter has numeric cm markings
- Range of sizes 8 Fr - 18 Fr
- Biocompatible, long lasting polyester umbilical tape
- Supplied with water soluble lubricant
- Colour coded clips for Fr size
- Single patient use
- Use for up to 4 weeks
- Supplied with spare clip and opening device

¹ Lynch A, et al. *Aus J Otolaryngol* 2018;1:8;

² CORGRIP* NG / NI feeding tube retention system, Instructions for Use, 2017.

CORFLO* PEG Kit

Each kit contains the following components:

- Polyurethane PEG tubing with traction removable internal bolster.
- Ring shaped internal retention bolster (see table for conical shaped version)
- 2 types of external retention discs
- Fast release clamp
- ENFit® feeding connector
- Towel and drape
- Scalpel
- Jelly
- Stainless steel puncture needle and trocar G18

The "E1-type" of kits have:

- An additional G14 puncture needle with teflon trocar
- No jelly

The "E2-type" of kits have:

- A blue coated pull wire
- 1 G18 puncture needle, except Fr 12 kit has a Gauge 14 puncture needle with teflon coated trocar
- No scalpel

The product item#

50-5020E2 contains in addition to E2-type of kits:

- 3 needles: Gauge 14, 18, 25
- Scalpel, syringe and gauzes



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO* PEG Kit 12Fr, traction removable	50-6012E1	Box of 2
CORFLO* PEG Kit 16Fr, traction removable	50-6016E1	Box of 2
CORFLO* PEG Kit 20Fr, traction removable	50-6020E1	Box of 2
CORFLO* PEG Kit 12Fr, traction removable	50-6012E2	Box of 2
CORFLO* PEG Kit 16Fr, traction removable	50-6016E2	Box of 2
CORFLO* PEG Kit 20Fr, traction removable	50-6020E2	Box of 2
CORFLO* PEG Kit 12Fr, traction removable	50-6512	Box of 2
CORFLO* PEG Kit 16Fr, traction removable	50-6516	Box of 2
CORFLO* PEG Kit 20Fr, traction removable	50-6520	Box of 2
CORFLO* PEG Kit 20Fr, traction removable, conical shaped bolster	50-5520	Box of 2
CORFLO* PEG Kit 20Fr, traction removable, conical shaped bolster	50-5020E2	Box of 2
CORFLO* PEG Kit Over the Wire 16Fr, traction removable	50-4016E1	Box of 2
CORFLO* PEG Kit Over the Wire 20Fr, traction removable	50-4020E1	Box of 2

CORFLO[®] PEG Jejunal Tube, polyurethane, radiopaque, non-sterile



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
Jejunal tube for use with CORFLO [®] PEG Fr 16: 6Fr x 91cm, 1 gr weighted, pull along technique	50-7366	Box of 5
Jejunal tube for use with CORFLO [®] PEG Fr 20: 8Fr x 91cm, 3 gr weighted, pull along technique	51-7368	Box of 5
Jejunal tube for use with CORFLO [®] PEG Fr 20: 10Fr x 91cm, 3 gr weighted, pull along technique	50-7361	Box of 5
Jejunal tube for use with CORFLO [®] PEG Fr 20: 8Fr x 109cm, unweighted, over-the-wire technique	50-1438	Box of 2
Jejunal tube for use with CORFLO [®] PEG Fr 20: 10Fr x 109cm, unweighted, over-the-wire technique	50-1431	Box of 2

CORFLO[®] PEG Adapter Repair Kit



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO [®] PEG Adapter Repair Kit 12Fr	50-6112	Box of 5
CORFLO [®] PEG Adapter Repair Kit 16Fr	50-6116	Box of 5
CORFLO [®] PEG Adapter Repair Kit 20Fr	50-6120	Box of 5

ENTERAL ACCESSORIES

CORFLO* PEG Accessories



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO* PEG Silicone Fixation Device 12Fr	SL-0012	Box of 10
CORFLO* PEG Silicone Fixation Device 16Fr	SL-0016	Box of 10
CORFLO* PEG Silicone Fixation Device 20Fr	SL-0020	Box of 10
CORFLO* PEG Retention Device 12Fr	RA-0012	Box of 10
CORFLO* PEG Retention Device 16Fr	RA-0016	Box of 10
CORFLO* PEG Retention Device 20Fr	RA-0020	Box of 10
CORFLO* PEG Quick Release Clamp for 12Fr & 16Fr	CL-1216	Box of 20
CORFLO* PEG Quick Release Clamp for 20Fr	CL-0020	Box of 20



CORFLO* Tube Accessory

PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORFLO* Extension Set	40-1007	Box of 30

FARRELL* Valve System



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
FARRELL* Valve System: Gastric Pressure Relief System - The FARRELL* Valve provides a channel to constantly decompress the stomach, allowing the stomach to fill at its own pace	44-4100	Box of 30

aspHirate pH Indicator Strips

Results from an audit examining the viability of pH indicator strips for nasogastric tube placement found participants made up to 11.15% errors reading the strips and reported anxiety in their use.¹

Study results supported the development of new pH indicator strips that are more usable and less likely to be misread.¹

pH test strips for the semi-quantitative indication of pH in human gastric tube aspirate for the assessment of correct tube placement

INNOVATING THE STANDARD - IMPROVING PATIENT SAFETY

- Easy to read at critical clinical decision values
- Accuracy supported with independent test data
- Adjusted to read water at pH 7
- Unique scale for accurate documentation
- Unique single scale printed lengthwise for ease of use
- Pictorial instructions alongside scale
- Unique 2 colour match
- Desiccant lid



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
aspHirate pH Indicator Strip: unique 2 colour pad and scale, pictorial IFU for ease of use, 0.5 increments, CE IVD for human gastric aspirate 3.5 - 7.5	code 150.1	100 per pot

Enteral Syringe, single use

ENFit® Enteral Feeding Syringes and Accessories from Avanos are designed to increase patient safety and optimal delivery of enteral feeding by reducing the risk of tubing misconnections. They are ISO 80369-3 compatible and do not contain natural rubber latex, BPA or DEHP.



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
Enteral Low-Dose Tip Syringe with ENFit® Compatible Connector, 1ml	SYR-01S	Box of 100
Enteral Low-Dose Tip Syringe with ENFit® Compatible Connector, 3ml	SYR-03S	Box of 100
Enteral Low-Dose Tip Syringe with ENFit® Connector, 5ml	SYR-05S	Box of 100
Enteral Syringe with ENFit® Connector, 10ml	SYR-10S	Box of 100
Enteral Syringe with ENFit® Connector, 20ml	SYR-20S	Box of 50
Enteral Syringe with ENFit® Connector, 60ml	SYR-60S	Box of 50

¹ Ref BorsciS, et al. *PLoS ONE* 2017;12(11):e0189013

CORSTOP* a.c.e. stopper

- CORSTOP* is indicated for patients in e.g. antegrade continence enemas (a.c.e.) or the urine bladder (Mitrofanoff procedure) where they are required to seal the tract.
- It is used to maintain patency in between insertions of e.g. enema catheter.
- Also indicated as gastrostomy stopper in case of accidental removal of a gastrostomy feeding tube to prevent closing or narrowing of the stoma tract.



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CORSTOP* a.c.e. stopper 8Fr x 15mm	CS8-15	Each
CORSTOP* a.c.e. stopper 8Fr x 30mm	CS8-30	Each
CORSTOP* a.c.e. stopper 8Fr x 60mm	CS8-60	Each
CORSTOP* a.c.e. stopper 10Fr x 15mm	CS10-15	Each
CORSTOP* a.c.e. stopper 10Fr x 30mm	CS10-30	Each
CORSTOP* a.c.e. stopper 10Fr x 60mm	CS10-60	Each
CORSTOP* a.c.e. stopper 12Fr x 30mm	CS12-30	Each
CORSTOP* a.c.e. stopper 12Fr x 60mm	CS12-60	Each
CORSTOP* a.c.e. stopper 12Fr x 100mm	CS12-100	Each
CORSTOP* a.c.e. stopper 14Fr x 30mm	CS14-30	Each
CORSTOP* a.c.e. stopper 14Fr x 60mm	CS14-60	Each
CORSTOP* a.c.e. stopper 14Fr x 100mm	CS14-100	Each
CORSTOP* a.c.e. stopper hydrogel dressing	CS-01	Box of 30

CLOG ZAPPER* Enteral Feeding Tube Declogging System



PRODUCT DESCRIPTION	PRODUCT CODE	PACKAGED
CLOG ZAPPER* Enteral Feeding Tube Declogging System, for the treatment of feeding tube occlusion	090120090	10 per case

For Healthcare professional use only.

**From the UK, please call Avanos
at 0800 917 65 85 - fax 0800 169 02 35
or send an email to
customerservice.uk.ie@avanos.com**

**From IRL, please call Uniphar Hospital Services
at +353 1 468 84 56 - fax +353 1 466 07 82
or send an email to unipharhospitalservices@uniphar.ie**

**From N-Ireland - please call Sangers Surgical
at 028 9056 7135 - fax 028 9070 5623
or send an email to Hospital@sangersaah.co.uk**

**From all other countries, please call Avanos
at +32 2 700 68 51 - fax +32 2 711 26 91
or send an email to customerservice.export@avanos.com**

AVANOS

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. The COLOR YELLOW is a Registered Trademark or Trademark of Avanos Medical, Inc. or its affiliates. © 2019 AVNS. All rights reserved.

ENFit® is a Registered Trademark or Trademark of Global Enteral Device Supplier Association, Inc. Used with their permission. COPY-05145